



1863/64 Fuse Part Numbers			
GenRad/IET Part Number	Current Rating	Case Size	Type
 <p style="text-align: right;">Older Style with two holders Fuses are item 4 and 5.</p>			
5330-0300-00	1/16 A (0.063A)	1 ¼ (1.25")	Slo-Blo
5330-0450-00	1/8 A (0.125A)	1 ¼ (1.25")	Slo-Blo
5330-0700-00	1/4 A (0.250A)	1 ¼ (1.25")	Slo-Blo
 <p style="text-align: right;">Newer Style with IEC Inlet Module</p>			
Recommend Use 0034.3111	1/16 A (0.063A)	5mm x 20mm	Slo-Blo
Recommend Use 0034.3111	1/8 A (0.125A)	5mm x 20mm	Slo-Blo
0034.3111 (For Use from 90V to 240VAC)	1/4 A (0.250A)	5mm x 20mm	Slo-Blo